Precious Metals Assaying Laboratories

Laboratorio analisi metalli preziosi

Analytical scale and anti-vibration table

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Assays bench with crushing, rolling and annealing system

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Separation Hood

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Fire assay method for PM's determination is an old well known analytical method which is clearly described in ISO 11426-1998.

For gold assays by cupellation, the alloys are inquarted with silver, compounded with lead (and Cu for high purity alloys) and cupelled in a cupellation furnace until a precious metal button is obtained. After flattening and rolling, silver is extracted in nitric acid and gold weighed.

Following a brief outline of the process and the different steps:

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Operation:	Requested equipment:
Sampling	Precision bench drill
Weighting	Analytical scale and anti vibration table
Inquartation	Assays bench
Cupellation	Cupellation Furnace + Cartridge Filter for Lead Oxide removal
Milling operation and Annealing	Assay bench
Chemical attack by nitric acid	Separation hood
Weighting and calculation	Analytical scale and anti vibration table

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SEPARATION HOODS 6 OR 12 PLACES

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Separation hoods 6 or 12 places

Machines used to separate the different metals by nitric acid reaction.

Made of polypropilene, equipped with arrays of electric heaters for the assaying with glass flasks.

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Analytical scale and anti vibration table

High precision analytical scale used to weigh the sample to be assale.

Provided with self calibration.

An anti-vibrating table with marble slab can also be provided on request.

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Furnaces

Cupellation furnace with heated resistor or by silite rods.

Extractable muffle made with silicon carbide (carborundum)

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Assays bench

Stainless steel work top with anvil to crush cupels ball. Available with manual or motorie rolling mill with foot pedal control.

Electric annealing system provided with voltage regulator for the gold foil treatment.

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Cartridge filter for lead oxide

With furnaces for cupellation we recommend to purchase a filter for lead oxides coming from cupellation. The cartridge cleaning is automatically done, controlled by a timer. The plant is equipped with a bag for recovery powder.

The plant is available with a bag for powder recovery.

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LABORATORY HOODS

Balestri Technologies designs and realizes lab fume hoods made in different materials (PP, PPs) complete of vacuum system and customized on customer's requirements.

The usage of materials and high quality components guarantees the maximum security, according to the actual code and regulations, and long lasting durability in time. Our fume hoods are very versatile; they can easily equip industrial, pharmaceutical, scientific or scholastic labs.

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